# Stieber Rippenbogen Retractor

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
## Teile

G-110 Stieber Kit

- A Optional: LEE Center
- B Optional: Hydra Arm
- C Vertical Pole with Wingnut Haltestab vertikal Gewinde
- D Vertical Pole with V-Notch Haltestab vertikal Gabel
- E Two-Piece Horizontal Bar Querstange verschraubt
- F Grooved 3/D Swinger Rail Clamp (x2) Socket
- G Stieber Small Platform (Left & Right) Plattform klein links & rechts
- H Stieber Large Platform (Left & Right) Plattform groß links & rechts
- I Stieber Turning Handle (x2) Kurbel 2 mal je links & rechts identisch
- J Stieber Needle Nadelöffnung
- K 2 Pairs of Stieber Ratcheted Handle Zahnstangen mit Retractoren 2x

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
### Assemble the Platforms


**1** Place the STIEBER TURNING HANDLE through the hole of the SMALL RIGHT PLATFORM.

Kurbel in kleine Plattform.

**2** Insert the end of a RATCHETED BLADE **WITH THE BLADE CURVING DOWN** into the SMALL RIGHT PLATFORM and feed through by circling the HANDLE clockwise.

Zahnstangen in kleine Plattform.

**3** Fit the hole in the STIEBER SMALL RIGHT PLATFORM over the pin on the STIEBER LARGE RIGHT PLATFORM. Fold the pin down to lock into place.

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### Install the Frame

**4** Fully open the GROOVED 3/D SWINGER RAIL CLAMP by turning the WINGNUT counterclockwise.

Socket komplett öffnen.


**5** Attach the GROOVED 3/D SWINGER RAIL CLAMP onto the bed rail on the **LEFT SIDE OF THE PATIENT**. Hook the larger and wider jaw onto the bed rail edge with the lower jaw resting underneath the bed rail.

Socket an OP Schiene.


**6** Insert the VERTICAL POLE WITH WINGNUT through the top of the GROOVED 3/D SWINGER RAIL CLAMP with the **POLE CURVING AWAY** from the patient. Vertikale Stange einsetzen. Biegung weg vom Patient,


**IMPORTANT**  
The VERTICAL POLE WITH WINGNUT **MUST BE SET UP ON THE PATIENT'S LEFT SIDE.**

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**7** Move the INTERLOCKING SLEEVE to sit on top of the GROOVED 3/D SWINGER RAIL CLAMP. Align the sleeve notch with the retaining pin to secure. Secure the GROOVED 3/D SWINGER RAIL CLAMP by tightening the WINGNUT clockwise.

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
**9** Fully open the GROOVED 3/D SWINGER RAIL CLAMP by turning the WINGNUT counterclockwise.

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
**10**

Attach the **GROOVED 3/D SWINGER RAIL CLAMP** onto the bed rail on the **RIGHT SIDE OF THE PATIENT**. Hook the larger and wider jaw onto the bed rail edge with the lower jaw resting underneath the bed rail. Sockel an OP Schiene.


**11**

Insert the **VERTICAL POLE WITH V-NOTCH** through the top of the **GROOVED 3/D SWINGER RAIL CLAMP** with the **POLE CURVING AWAY** from the patient. Vertikale Stange einsetzen. Biegung weg vom Patient,

**IMPORTANT**  
MAKE SURE THE  
POLE IS CURVING  
AWAY FROM THE  
PATIENT.

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**12**


Move the **INTERLOCKING SLEEVE** to sit on top of the **GROOVED 3/D SWINGER RAIL CLAMP**. Align the sleeve notch with the retaining pin to secure. Secure the **GROOVED 3/D SWINGER RAIL CLAMP** by tightening the **WINGNUT** clockwise.

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
**17**  
Set the end of the HORIZONTAL BAR onto the VERTICAL POLE WITH V-NOTCH.

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**18**  
Attach the HYDRA YOKE DOUBLE ARM to the LEES CENTER PLATFORM making sure **THE SCALLOPED EDGE OF THE YOKE FACES DOWN**, and secure with the WINGNUT.  
OPTIONAL: Hydra Arm aufsetzen.


**19**  
Secure the HORIZONTAL BAR to the VERTICAL POLE WITH WINGNUT with the set aside WINGNUT.

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### Configure the Arms

**20**  
Insert a BLADE HANDLE into the WRIST JOINT and secure by turning the WRIST WINGNUT clockwise.  
OPTIONAL: Hydra mit Retractoren bestücken.

**21**  
Position the HYDRA YOKE ARMS as desired.

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
**22**  
Secure the HYDRA ARMS by tightening the ELBOW WINGNUTS.

Loosen the ELBOW WINGNUTS to allow manipulation of the HYDRA ARMS.

OPTIONAL: Hydra fixieren.

### Stieber Needle

Use the STIEBER NEEDLE to assist in tightening LEES CENTER PLATFORM and STIEBER LARGE PLATFORMS if desired.

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